**Best. Nr.** CK 121  
**Inhalt:** 20 Tests

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HK = Hexokinase, G6P-DH = Glucose-6-Phosphat-dehydrogenase, ADP = Adenosin-diphosphat, ATP = Adenosin-triphosphat

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**Safety information**  
 The buffer solution (round cuvette) contain  
 - Sodium azide (<0.1 %). Do not swallow and avoid contact  
 with skin and mucous membranes.  
 - 0.6% Imidazole and is categorized as a dangerous  
 preparation according to EC Directives:  
 H360D: May damage the unborn child  
 Observe the safety advice on the packaging.  
 A safety data sheet is available on request.<sup>3)</sup>

**Storage and shelf life**  
 The reagents have to be kept at a temperature between  
 +2°C and +8°C until the expiry date indicated on the  
 packaging.  
 Stability of ready-to-use reagent:  
 at +2°C to +8°C 48 h  
 at +20°C to +25°C 8 h

**Measurement conditions**  
**Measurement device:** Vario Photometer II Diaglobal  
**Meas. wavelength:** 365 nm  
**Temperature:** 37°C  
 Additionally needed:  
 Mini centrifuge, Dry block heater (37°C)

**Measurement range**  
 Up to 1000 U/L  
 In case of exceeding this value, dilute the sample 1 + 4  
 with physiological saline solution, multiply the result by 5.

**Reference values**  
 Male < 170 U/L  
 Female < 145 U/L

**Preparation**  
 The dry block heater should be activated at least 20 min.  
 before measurement.

**Working instructions**

Pipette into round cuvettes:	
	Analysis
Enzyme substrate solution	500 µL
Pipette into ready-to-use reagent:	
Serum / plasma	20 µL
Mix thoroughly and put immediately into dry block heater.	

Up to six samples can be measured at the same time.  
 • Pipette 500µL enzyme substrate solution in cuvette  
 • with buffer solution (=ready-to-use reagent)  
 • Select <CK 121>  
 • Pipette 20µL serum / plasma in cuvette with ready-to-  
 use reagent, mix well  
 • Insert cuvette into 37°C dry block heater and start  
 measurement at the photometer  
 • Follow the operator guidance showing on display  
 • After finishing the measurements, read results in (U/L)  
 with arrow keys  
 • Remove cuvette

**Calculation**  
 Activity (U/L) of creatine kinase in plasma will be calculated.  
 $CK (U/L) = F \times \Delta A$   
 The computation formula is stored in the Photometer.


**Quality assurance**  
 For quality assurance we recommend universal control  
 sera from company Roche, www.roche.de:  
 PreciControl ClinChem Multi 1 / Multi 2 (4 x 5 mL)  
 Order-No.: 05 947 626 190 / 05 947 774 190  
 Ref.: Roche / Hitachi analyzers, Method: IFCC liquid

**Advice**  
 For CK measurement in capillary blood, we recommend  
 our test CK 321.

**Summary**  
**Biochemical fundamentals**  
 Creatine Kinase has a key role in the catabolism of muscle  
 cells and is responsible for reduction of ADP to ATP.  
 Measurable total activity of CK in serum consists of activity  
 of iso enzyme CK-MB, CK-MM and CK-BB. In blood plasma  
 of a healthy person a very low concentration of CK is normally  
 found. At muscle deterioration CK transgresses into blood  
 and blood CK level increases. Since evacuation from cell  
 occur by lymph channel the concentration climax will be  
 reached mostly after 6 - 8 hours.

**Training and competition control**  
 Increased CK level due to physical activity is an expression  
 of overstraining or unaccustomed strain of individual muscle  
 groups.  
 If the CK value increases to 900 U/L (15 µmol/Ls) during  
 training, a decrease of training exercise is recommended or  
 else overtraining is risked.<sup>4)</sup>  
 After extreme physical stress (marathon, triathlon) CK value  
 can increase beyond 4000 U/L (70 µmol/Ls).<sup>4)</sup>

**Measurement principle**  
 With Diaglobal test CK 121 total CK is measured.  
 Determination is based on a 37°C standard method of IFCC  
 (International Federation of Clinical Chemistry).<sup>5)</sup>

HK = Hexokinase, G6P-DH = Glucose-6-Phosphate-dehydrogenase,  
 ADP = Adenosine-diphosphate, ATP = Adenosine-triphosphate

The increase of NADPH is directly proportional to the CK  
 activity and is measured at 365 nm.

Note: Older technical literature often gives values that were determined  
 with the help of the older 25°C method. Values, measured with the new  
 IFCC method, average higher about a factor of 1.6.

**Performance parameters**  
**Specificity / interferences<sup>6)</sup>**  
 Falsely increased values due to strong haemolyse (Hb >2.0  
 g/L), no significant interferences due to bilirubin (up to 15  
 mg/dL).  
 Strong lipemic samples can cause problems due to high  
 absorptions. IN this case dilute the sample 1 + 4 with  
 physiological saline solution.

**Inaccuracy**  
 The reproducibility was checked using human and control  
 samples.

In series [n = 20]	Average [U/L]	Standard deviation [U/L]	CV [%]
Sample 1 Sample 2	156 502	4.84 9.54	3.1 1,9
From day to day [n = 20]	Average [U/L]	Standard deviation [U/L]	CV [%]
Sample 1 Sample 2	152 499	5.32 12.0	3.5 2.4

**Analytic sensitiveness**  
 Lower detection limit: 20 U/L (0.33 µkat/L)

**Comparison of methods**  
 A comparison of the Diaglobal test CK 121 (y) and another  
 commercial available test (x) resulted in the following  
 correlation according to the linear regression:  
 $y = 0.999x + 0,8$   
 $r = 0.994$   
 $n = 40$   
 Concentration range: 40 - 2000 U/L

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