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
HK = Hexokinase, G6P-DH = Glucose-6-Phosphat-dehydrogenase,
 ADP = Adenosin-diphosphat, ATP = Adenosin-triphosphat

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Order no. CK 321
Contents: 20 tests
 20 pipette tips, 500µL

Method
 IFCC 37°C^{1,2)} / DGKC modified

Sample material
 Capillary or EDTA blood.
 Set capillary blood immediately in reaction tube "R".
 Stability of the sample in reaction tube "R":
 6 hours (+15°C to +25°C)

Reagents
 Contents / concentrations:
 1. Starter reagent (caps in PE-bottle)
 ADP 2 mmol/L, G6-PDH >2 KU/L
 2. Buffer solution (pre-portioned in round cuvettes)
 Imidazole buffer 100 mmol/L, pH = 6.00 (37°C), Glucose 20 mmol/L,
 Magnesium acetate 10 mmol/L, EDTA 2 mmol/L, AMP 5 mmol/L,
 NADP 2 mmol/L, N-Acetyl-cysteine 20 mmol/L, Diadenosin-penta-
 phosphate 10 µmol/L, HK > 4.0 KU/L; Sodium azide < 0.1 %
 3. Solution for dilution of blood (reaction tubes "R")
 Creatine phosphate 30 mmol/L, Sodium chloride 3 g/L,
 Sodium azide < 0.1 %

Safety information
 The buffer solution (round cuvette) contain
 - Sodium azide (<0.1 %). Do not swallow and avoid contact
 with skin and mucous membranes.
 - 0.6% Imidazole and is categorized as a dangerous
 preparation according to EC Directives:
 H360D: May damage the unborn child
 Observe the safety advice on the packaging.
 A safety data sheet is available on request.³⁾

Storage and shelf life
 The reagent have to be kept in a dark place at a
 temperature between + 2°C to + 8°C until the expiry date
 indicated on the packaging.

Measurement conditions
 Measurement device: Vario Photometer II Diaglobal

Meas. wavelength: 365 nm

Temperature: 37°C

Additionally required:
 Mini centrifuge, Dry block heater (37°C)

Measurement range
 Up to 1000 U/L
 In case of exceeding these value, repeat determination with
 sample volume of 20 µL (instead of 60 µL). Multiply the
 result by 3

Reference values
 Male < 170 U/L
 Female < 145 U/L

Preparation
 Reagents are ready to use.
 The dry block heater should be activated at least 20 min.
 before measurement.

Working instructions
 Up to six samples can be measured at the same time.
 • Withdraw 60µL capillary blood from finger pad or earlobe
 with end-to-end capillary and insert into reaction tube "R"
 • Transfer blood into solution by mixing strongly
 • Insert reaction tube "R" into Minicentrifuge and centrifuge
 for 1 minute
 • Pipette 500µL supernatant into round cuvettes (blanks)
 • Select <CK 321>
 • Screw caps from PE-bottle onto prepared round cuvettes
 and mix very strongly to solve the starter reagent
 • Insert cuvettes into 37°C dry block heater and start
 measurement at the same time
 • Follow the operator guidance showing on display
 • Upon completion of measurement: Adjust HCT value with
 arrow keys and confirm with <ENTER> *)
 • After finishing the measurements, read results in U/L with
 arrow keys
 • Remove cuvette

*) HCT value of 40 % is set as standard value.
 Please note that in case of measuring CK 321 with control
 serum the HCT-value has to be keyed in as 0%.

Calculation
 Activity (U/L) of Creatine kinase in plasma will be
 calculated. For exact measurement results the haematocrit
 values has to be considered.

$$CK (U/L) = F \times \Delta E / (1 - 0.01 \times HCT)$$

F = calibration constant, HCT = Haematocrit in %
 The algorithm for calculating the results is stored in the
 Vario Photometer II.

Conversion: µkat/L x 60 = U/L

Quality assurance
 For quality assurance we recommend universal control sera
 from company Roche, www.roche.de:
 PreciControl ClinChem Multi 1 / Multi 2 (4 x 5 mL)
 Order-No.: 05 947 626 190 / 05 947 774 190
 Ref.: Roche / Hitachi analyzers, Method: IFCC liquid

Summary
Biochemical fundamentals
 Creatine Kinase has a key role in the catabolism of muscle
 cells and is responsible for reduction of ADP to ATP.
 Measurable total activity of CK in serum consists of activity
 of iso enzyme CK-MB, CK-MM and CK-BB. In blood plasma
 of a healthy person a very low concentration of CK is normally
 found. At muscle deterioration CK transgresses into blood
 and blood CK level increases. Since evacuation from cell
 occur by lymph channel the concentration climax will be
 reached mostly after 6 - 8 hours.

Training and competition control
 Increased CK level due to physical activity is an expression
 of overstraining or unaccustomed strain of individual muscle
 groups.
 If the CK value increases to 900 U/L (15 µmol/Ls) during
 training, a decrease of training exercise is recommended or
 else overtraining is risked.⁴⁾
 After extreme physical stress (marathon, triathlon) CK value
 can increase beyond 4000 U/L (70 µmol/Ls).⁴⁾

Measurement principle
 With Diaglobal test CK 321 total CK is measured.
 Determination is based on a 37°C standard method of IFCC⁵⁾
 (International Federation of Clinical Chemistry).

HK = Hexokinase, G6P-DH = Glucose-6-Phosphate-dehydrogenase,
 ADP = Adenosine-diphosphate, ATP = Adenosine-triphosphate

The increase of NADPH is directly proportional to the CK
 activity and is measured at 365 nm.

Note: Earlier technical literature often gives values that were
 determined with the help of the 25°C method. Values, measured with
 the new IFCC method, average higher about a factor of 1.6.

Performance parameters
Specificity / interferences⁶⁾
 Falsely increased values due to strong haemolysis (Hb>2.0
 g/L), no significant interferences due to bilirubin (up to 15
 mg/dL).
 Strong lipemic samples can cause problems due to high
 absorptions.

Inaccuracy
 The reproducibility was checked using human and control
 samples.

In series [n = 20]	Average [U/L]	Standard deviation [U/L]	CV [%]
Sample 1	150	4.96	3.3
Sample 2	503	12.3	2.4
From day to day [n = 20]	Average [U/L]	Standard deviation [U/L]	CV [%]
Sample 1	148	5.48	3.7
Sample 2	508	14.2	2.8

Analytic sensitiveness
 Lower detection limit: 20 U/L (0.33 µkat/L)

Comparison of methods
 A comparison of the Diaglobal test CK 321 (y) and another
 commercial available test (x) resulted in the following
 correlation according to the linear regression:
 $y = 0.984x - 3.65$
 $r = 0.996$
 $n = 45$
 Concentration range: 40 - 2600 U/L

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