Best. Nr. GOT 442
Inhalt: 40 Tests

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Method
 UV test, IFCC 37°C^{1,2}, modified,
 without pyridoxal phosphate activation,
 working process with sample start

Sample material
 Serum, heparin or EDTA plasma
 activity diminution in serum
 at +2°C to +8°C: after 3 days ca. 10 %
 at +20°C to +25°C: after 3 days ca. 17 %
 Haemolysis disturbs.

Reagents
 Content / concentrations:
 1. Starter reagent (caps in PE-bottle)
 NADH > 0.18 mmol/L, α-ketoglutarate 13.3 mmol/L
 2. Enzyme substrate solution
 (pre-portioned in round cuvettes)
 Lactate dehydrogenase (LDH) from pig muscle > 0.9 kU/L, Malate
 dehydrogenase (MDH) from pig heart > 0.6 kU/L,
 L-aspartate 240 mmol/L, Sodium azide < 0.1 %, TRIS buffer
 pH 7.65 (37°C) 80 mmol/L

Safety information
 The buffer solution (round cuvette) contains sodium azide
 (<0.1 %). Do not swallow and avoid contact with skin and
 mucous membranes. If desired a safety data sheet will be
 provided.³⁾

Storage and shelf life
 Reagents can be kept at a temperature between +2°C and
 +8°C until the expiry date indicated on the packaging.

Measurement conditions
 Measurement device: Vario Photmeter II
 Meas. wavelength: 365nm
 Temperature: 37°C

Measurement range
 10 - 500 U/L
 Should values exceed this range, dilute serum/plasma 1+10
 with physiological NaCl solution. Pipette 0.1 mL serum to
 1.0 mL NaCl solution, mix and repeat the determination with
 50 µL of this dilution. Result x 11.

High activity sera may show a low starting absorption if the
 NADH has already been consumed for the greatest part at
 the beginning of the reaction. (Display: ABS 1 too low). In
 this case, dilute sample as indicated above.

Preparation
 Switch on the dry block heater at least 20 minutes before
 measurement.

Working instructions
 You may measure up to six samples simultaneously.
 • Incubate GOT 442 cuvettes 7 minutes at 37°C.
 • Pipette 50 µL serum/plasma with end-to-end capillary into
 cuvettes GOT 442, *not mixing yet!*
 • Screw caps from the PE bottle onto the cuvettes and
 dissolve the caps' content by inverting several times, mix
 thoroughly. Hereby the sample has to be streamed out of
 the capillary completely.
 • Place back cuvettes into the 37°C thermostat.
 • Select test <GOT>, start measurement.
 • Follow user instructions on display.
 • The results (U/L) are shown on display by pressing the
 right arrow key after the measurements.

Calculation
 The activity (U/L) of GOT in plasma becomes calculated.

$$\text{GOT (37°C)} = F \times \Delta\text{ABS} - 3.2$$

F = Calibration constant
 The calculation formula is coded in the Vario Photometer II.

Conversion in µkat/L: U/L x 0.0167 = µkat/L
 Conversion in U/L (25°C):
 U/L (37°C) x 0.48 = U/L (25°C)

Quality assurance
 For quality assurance we recommend universal control sera
 from company Roche, www.roche.de:
 PreciControl ClinChem Multi 1 / Multi 2 (4 x 5 mL)
 Order-No.: 05 947 626 190 / 05 947 774 190
 Ref.: Roche / Hitachi analyzers,
 Method: IFCC without pyridoxal phosphate

Reference values⁴⁾ for serum / plasma

	U/L 37°C	µkat/L 37°C
Men	until 35	until 0.59
Women	until 31	until 0.52

Summary^{1,2)}
 The glutamate oxalacetate transaminase, GOT (new
 nomenclature: aspartyl aminotransferase, AST) belongs to
 the transaminases, a group of enzymes abundantly found in
 particular in liver, myocardium and skeletal muscle and
 transferred into blood in case of any damage of these organs.

Indications / diagnostic significance:
 diagnostics and follow-up of
 - liver and biliary tract diseases
 - skeletal muscle diseases
 - heart attack

GOT and GPT are sensitive indicators for liver diseases.
 After a myocardial infarct the GOT activity in serum increases
 and reaches its highest value after 2 days. In Germany the
 optimised standard method of the German Society for
 Clinical Chemistry (DGKC) has been applied so far for GOT
 determination. By now this method has been replaced by the
 IFCC method 37°C.

Measuring principle^{2,4)}
 Under catalysis with GOT α-ketoglutarate becomes
 converted with L-aspartate into L-glutamate and oxalacetate.
 The following reduction of oxalacetate to L-malate, catalysed
 by MDH, causes NADH to oxidise to NAD⁺.

The NADH consumption is measured at 365 nm. The
 absorption decrease is proportional to the GOT activity in
 serum / plasma.

Performance parameters
Specificity / interferences
 Haemolysis disturbs as GOT is set free out of erythrocytes.
 No significant interference because of bilirubin and
 lipaemia. A very intense lipaemia may cause a too high
 starting extinction. In this case dilute sample 1 + 1 or 1 + 4.
 The determination is conducted without additional pyridoxal
 phosphate. Thus GOT values of samples with vitamin B6
 deficiency may result too low.

Inaccuracy
 The reproducibility was checked using human and control
 samples.

In series [n = 18]	Average [U/L]	Standard deviation [U/L]	VK [%]
Sample 1 Sample 2	42.2 146	1.8 3.6	4.3 2.5
From day to day [n = 18]	Average [U/L]	Standard deviation [U/L]	VK [%]
Sample 1 Sample 2	42.3 149	2.2 6.5	5.1 4.3

Analytic sensitiveness
 Lower detection limit: 10 U/L (0.16 µkat/L)

Comparison of methods
 A comparison of the Diaglobal test GOT 442 (x) and a
 commercially available test (y) resulted in the following
 correlation according to the Passing/Bablok⁵⁾ process:
 $y = 1.013x + 2.285$
 $r = 0.997$

n = 24
 Concentration range: 15 - 500 U/L

Bibliography
 1. Thomas L. Labor und Diagnose. 4th edition. Marburg: Die
 Medizinische Verlagsgesellschaft 1995: 120
 2. International Federation of Clinical Chemistry, Scientific
 Committee. Eur J Clin Chem Clin Biochem 1977; 15: 39-51
 3. http://www.diaglobal.de/de/service/downloads/index.html
 4. Schumann G et al. Clin Chem Lab Med 2002 ; 40 :725-733
 5. Passing H, Bablok W. A new biometric procedure for testing the
 equality of measurements from two different analytical methods.
 J Clin Chem Clin Biochem. 1983; 21:709-720