**Best. Nr. MG 013**

Inhalt: 2 x 50 mL Farbreagenz

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 Calcium-Magnesium-Standard Best. Nr. CA ST  
 Diaglobal-Rundküvetten, 40 Stk. Best. Nr. LH 075

**Methode**

 Xylidylblau<sup>1,2)</sup>
**Probenmaterial**

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Contents: 2 x 50 mL Colour reagent

Additionally needed:

 Calcium-Magnesium-Standard Order No. CA ST  
 Diaglobal-Round Cuvettes, 40 pcs. Order No. LH 075

**Method**

 Xylidyl blue<sup>1) 2)</sup>
**Sample material**

Serum, heparinised plasma or urine

 Do not use EDTA-plasma, haemolytical plasma or sera.  
 Add some drops of HCl to urine, dilute 1+4 with dest. Water.

Stability in serum: at +2°C to +25°C: 7 days

**Reagent**

Contents / concentrations of the colour reagent:

Xylidylblue 0.1 mmol/L; GEDTA 40 µmol/L; Potassium carbonate 70 mmol/L; Detergent 1.5%; Tris-buffer 200 mmol/L; Sodium azide &lt; 0.1%; Triton X-100 &lt; 1%, pH 11,0

**Safety information**

The reagent contains sodium azide (&lt; 0.1 %) and Triton X-100. Do not swallow and avoid contact with skin and mucous membranes.

 A safety data sheet is available on our website.<sup>3)</sup>
**Storage and shelf life**

The colour reagent can be kept at a temperature between +2°C and +8°C until the expiry date indicated on the packaging.

**Measurement conditions**

Measurement device: VET Photometer Diaglobal

Meas. wavelength: 520 nm

Temperature: Room temperature

**Measurement range**

0.10 - 2.10 mmol/L (0.24 - 5.11 mg/dL)

Should values exceed this range, dilute sample 1+1 with physiological saline solution.

Multiply the result by 2.

**Notes**

1. Protect reagent from sunlight. Close the bottle after using.
2. The magnesium assay is very sensitive. We recommend the application of disposables.
3. For the "on-site" preparation of plasma we recommend the Diaglobal-Mini centrifuge (Order No. DZ 002)

**Working instructions**

Pipette into round cuvettes:			
	Blank	Standard	Analysis
Colour reagent	1500 µL	1500 µL	1500 µL
Standard	--	10 µL	--
Sample	--	--	10 µL

Mix well and measure after 5 minutes at room temperature standard and analysis against blank.

- Select the <MG> test
- Insert cuvette with reagent into photometer (blank value)
- Remove blank cuvette after signal ton
- Insert cuvette with standard into photometer (standard-absorption is shown on the display)
- Remove cuvette
- Insert cuvette with sample (Analysis) into photometer
- Read the result

**Tip**

It is possible to prepare several vials for measurement at the same time by pre-pipetting 1500 µL of colour reagent into the appropriate number of Diaglobal round cuvettes.

Shelf life of the pipetted reagents in the vial: 5 hours

**Calculation**

Concentration of magnesium in serum/plasma/urine:

 $c \text{ [mmol/L]} = 0.8226 \times \text{Abs. analysis} / \text{Abs. standard}$ 
 $c \text{ [mg/dL]} = 2.0 \times \text{Abs. analysis} / \text{Abs. standard}$ 

Conversion in mg/dL: mmol/L x 2.4312 = mg/dL

**Quality assurance**

For quality assurance we recommend our calcium - /magnesium - control solution CA QS.

**Reference values for cattle<sup>4)</sup>**

	mmol/L	mg/dL
In serum / plasma	0.8 - 1.3	1.9 - 3.2

**Summary<sup>4)</sup>**

Magnesium is a component of many different enzymes which are involved in energy metabolism in the muscle cells. A reduced concentration of magnesium in the blood (hypomagnesaemia) leads to tetanic seizures, while an increase (hypermagnesaemia) leads to flaccid paralysis.

**Evaluation<sup>4)</sup>**

 Hypermagnesaemia
 

- Parturient paresis
- Renal insufficiency
- Exsiccosis
- Latrogenic oversupply

Hypomagnesaemia

- Chronic bowel diseases (malabsorption)
- Low-magnesium diet
- Acute pancreatitis
- Cholestasis
- Renal insufficiency
- Hypercalcaemia
- Hypocalcaemic ketoacidosis
- Hyperthyroidism
- Tetany (cow, sheep)

**Measurement principle<sup>1) 2)</sup>**

In alkaline solution with xylidyl blue, magnesium forms a colour complex which is measured at 520 nm. Calcium ions are masked with Glycol Ethylenediaminetetraacetic Acid (GEDTA).

The assessment is carried out using a standard.

 $\text{Mg}^{2+} + \text{Xylidyl blue} \rightarrow \text{Magnesium colour complex}$ 
**Performance parameters**
**Specificity / interferences**

No interference caused by calcium ions.

Bilirubin up to 20 mg/dL and lipidaemia do not interfere.

Do not use haemolytic serums as erythrocytes have high magnesium concentrations.

**Inaccuracy**

The reproducibility was checked using human and control samples.

In series [n = 20]	Average [mmol/L]	Standard deviation [mmol/L]	VK [%]
Sample 1	0.88	0.02	2.3
Sample 2	1.15	0.02	2.1
Sample 3	1.45	0.02	1.6
From day to day [n = 20]	Average [mmol/L]	Standard-deviation [mmol/L]	VK [%]
Sample 1	0.67	0.04	5.8
Sample 2	0.93	0.05	5.4
Sample 3	1.38	0.05	3.9

**Analytic sensitiveness**

Lower detection limit: 0.08 mmol/L (0.2 mg/dL)

**Comparison of methods**

A comparison of the Diaglobal test MG 013 (y) and a commercially available test (x) resulted in the following correlation according to the linear regression:

$$y = 0,982x + 0,02$$

$$r = 0,993$$

n = 26, concentration range: 0.5 - 2.0 mmol/L

**Information on disposal**

Waste code number 180106:

Vials with reagent are considered hazardous waste. Do not allow reagent to reach surface water or sewage system.

Dispose of in accordance with official regulations.

Non-contaminated and completely empty packaging can be recycled.

Non-contaminated and completely empty packaging can be recycled.

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