Best. Nr. TRI 742
Inhalt: 40 Tests

Methode
 Enzymatischer Farbstest, GPO-PAP-Methode

Probenmaterial
 Biodiesel

- Reagenz**
 Inhalt / Konzentrationen:
- Startreagenz (Kappen in PE-Flasche)
 L-Glycerin-3-phosphat-oxidase (GPO) aus Mikroorganismen > 3,5 kU/L, Glycerokinase (GK) aus *Bacillus stearothermophilus* > 0,9 kU/L, Peroxidase (POD) > 3,5 kU/L, ATP 2,4 mmol/L, 4-Aminophenazon 0,15 mmol/L
 - Pufferlösung (vorportioniert in Rundküvetten)
 Lipoproteinlipase aus Mikroorg. > 7,5 kU/L, 2,4-Dichlorphenol 4 mmol/L, Natriumazid < 0,1 %, Triton X-100 < 1 %, PIPES-Puffer 50 mmol/L, pH 7,5

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
Messprinzip

GK = Glycerokinase, GPO = Glycerin-3-phosphat-oxidase

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Literatur
 1. <http://www.diaglobal.de>

Order No. TRI 742
Contents: 40 tests

Method
 Enzymatic colorimetric test, GPO-PAP method

Sample material
 Biodiesel

Reagents
 Contents / concentrations of the ready-to-use solution:
 1. Starter reagent (screw caps)
 L-glycerol-3-phosphate-oxidase (GPO) from microorganisms
 > 3.5 kU/L, Glycerokinase (GK) from bacillus stearothermophilus
 > 0.9 kU/L, Peroxidase (POD) > 3.5 kU/L, ATP 2.4 mmol/L,
 4-Aminophenazone 0.15 mmol/L
 2. Buffer solution (pre-portioned in round cuvettes)
 Lipoprotein lipase from microorg. > 7.5 kU/L, 2,4-Dichlorophenol
 4 mmol/L, Sodium azide < 0.1 %, Triton X-100 < 1%, PIPES-buffer
 50 mmol/L, pH 7.5

Safety information
 The buffer solution (round cuvette) contains sodium azide
 (< 0.1 %) and Triton X-100. Do not swallow and avoid contact
 with skin and mucous membranes. A safety data sheet is
 available on our website.¹⁾

Storage and shelf life
 The test reagents can be kept at a temperature between
 +2°C and +8°C until the expiry date indicated on the
 packaging. Please take the screw caps out of the container
 just before the analysis.

Measurement conditions
 Measurement device: Biodiesel Photometer Diaglobal

Meas. wavelength: 520nm

Temperature: Room temperature

The algorithm to compute the analysis result is coded in the
 above-named photometers.

Measurement range
 Mode <TRI> 0.1 - 23.0 g/dL
 Mode <TRI conc> 1.1 - 165 g/dL

Working instructions
 The measurement can be performed as a single or serial
 measurement (with a balancing of the A(0) counts = blank
 values).

Single measurement

Pipette into round cuvette:	
	Analysis
Sample	1 µL
Add 1 µL sample with capillary pipette, consider ring mark: 1-2-3-4-5 µL Mix vertical and intensive for 30 seconds or mix with vibration mixer. Wait for 1 minute before use. Solution needs to be clear.	

- Select the <TRI> test
- Insert analysis cuvette (blank value)
- Screw the cap from PE-bottle onto the cuvette, dissolve
 the starter reagent by inverting several times
- Press [ON/ENTER]
- Insert cuvette into the photometer immediately
- Wait for result

Calculation
 Absorbance x factor = TRI [g/dL]

Dilute when exceeding the measurement range
 If the measurement range is exceeded (display: > 15 g/dL)
 dilute the sample with freed biodiesel 1+10 and repeat the
 measurement in mode <TRI conc> with 1 µL diluted sample.
See Biodiesel dilution procedure
 The displayed measurement value has already been
 multiplied by the factor 11 and may be read directly.

Tip
 We recommend our capillaries 1 - 5 µL for pipetting 1µL
 sample order no. LH 021 and the micropipettor (pipettor
 help) order no. LH 007

Measurement principle

Triglycerides + H ₂ O	Lipase	⇌	Glycerol + Free fatty acids
Glycerol + ATP	GK	→	Glycerol-3-phosphate + ADP
Glycerol-3-phosphate + O ₂	GPO	→	Dihydroxyacetone Phosphate + H ₂ O ₂
H ₂ O ₂ + 2,4-Dichlorophenol + 4-Aminophenazone	Peroxidase	→	Quinonimine dye

GK = Glycerokinase, GPO = Glycerol-3-phosphate-oxidase

The concentration of the quinonimine dye is a measure for
 the triglyceride concentration in blood and serum/ plasma
 respectively and is measured photometrically at 520 nm.
 Reaching the end point of the reaction is identified
 automatically by the device itself.
 Free glycerol in serum, whose concentration corresponds to
 a triglyceride value of 10 mg/dl if the person is healthy, is
 measured, too.

Calibration function
 Besides triglycerides also monoglycerides and diglycerides
 as well as free glycerine become measured.

In case of lower triglyceride concentration we recommend to
 adjust the result through a calibration function that has to be
 determined empirically. For this the unadjusted values,
 measured in g/dL, are plotted against values that have been
 measured parallel with the reference method (gas
 chromatography).

Calculation of total glycerine
 The unadjusted triglyceride value allows to calculate total
 glycerine (the sum of free and bound glycerine n. DIN EN
 14214) as well:

TRI (g/dL) / 877 x 92 = TotGly (g/dL)
 TotGly (g/dL) / 0.8776 = TotGly (Ma%)
 877 = Molecular weight triglycerides [g/mol]
 92 = Molecular weight glycerine [g/mol]
 0.8776 = Density biodiesel [g/cm³]

Biodiesel dilution procedure

- Pipette 500µL biodiesel into a mixing vessel
- Withdraw 50µL biodiesel sample with the capillary and
 insert into the mixing vessel
- Close mixing vessel and mix thoroughly
- Select mode <TRI conc> and realise the measurement as
 usual with 1 µL of the diluted sample from the mixing
 vessel (see instruction „Step by Step“)
- The displayed measurement value has already been
 multiplied by the factor 11.

Accessories for dilution of biodiesel
 500µL pipette with 50 tips, order no. LH 500
 Mixing vessels, 500 pcs., order no. LH 050
 50µL capillaries, 100 pcs., order no. LH 056

Information on disposal
 Waste code number 180106:
 Vials with reagent are considered hazardous waste. Do not
 allow reagent to reach surface water or sewage system.
 Dispose of in accordance with official regulations.
 Non-contaminated and completely empty packaging can be
 recycled.

Bibliography
 1. <http://www.diaglobal.de>