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Adaption instructions for the quantitative determination of free fatty acids in serum/plasma using the Vet Photometer

Order no. NEFA 013

Contents R1 : 4 x 50 mL*
R1a: 4 x for 50 mL*
R2 : 4 x 25 mL*
R2a: 4 x for 25 mL*

Additionally required:

NEFA standard (2x10 mL)* Order no. NEFA ST

Diaglobal round cuvettes, 40 pcs. Order no. LH 075

* Product from FUJIFILM Wako Chemicals Europe GmbH

Method

ACS-ACOD method

Testing principle

See FUJIFILM package insert.

Sample material

Human serum / plasma as well as serum/plasma from cattle and other large animals.

Stability in serum: 2 days (+4 °C). It is advisable to determine the values immediately after blood is taken, wherever possible, as various serum enzymes break down fats to form free fatty acids. Freeze samples that are to be stored.

Reagents

Contents/concentrations of the ready-to-use solution: see FUJIFILM package insert.

Safety notice

Please observe warnings and precautionary measures stated on the box or package insert.

Storage and shelf life

At +2 °C to +8 °C, the reagents will keep until the expiry date stated on the packaging.

Measurement conditions

Instrument: Diaglobal Vet Photometer

Measurement wavelength: 520 nm

Temperature: Room temperature

Note

For collecting plasma on-site, we recommend the Diaglobal minicentrifuge (order no. DZ 002).

Preparation of reagents

Reagent R1: Mix contents of bottle R1a with solution R1 (50 mL).

Reagent R2: Mix contents of bottle R2a with solution R2 (25 mL).

Store the ready-to-use reagents between +2 °C and +8 °C and use within one month.

Working instructions

Pipette into round cuvettes:			
	Blank	Standard	Analysis
Reagent R1	1000 µL	1000 µL	1000 µL
Standard	--	50 µL	--
Sample	--	--	50 µL
Mix and leave to stand for 10 minutes at RT.			
Reagent R2	500 µL	500 µL	500 µL
Mix and leave to stand for 10 minutes at RT. Measure "Standard" and "Analysis" against "Blank".			

Procedure

- Select the <NEFA> test.
- Insert "Blank" cuvette into photometer (zero-point adjustment).
- Remove cuvette after the beep.
- Insert "Standard" cuvette into photometer (the display shows the standard absorbance).
- Remove cuvette.
- Insert cuvette containing sample (= "Analysis") into photometer; the measurement result is displayed in mmol/L.

Calculation

The algorithm for the calculation is stored in the device.

NEFA concentration c in serum/plasma:

$c [\text{mmol/L}] = 1.00 \times \text{Abs. Analysis} / \text{Abs. Standard}$

Conversion to mg/dL: $\text{mmol/L} \times 28.25 = \text{mg/dL}$

Test performance characteristics

See FUJIFILM package insert.

Information on disposal

Waste code number 180106:

Vials with reagent are considered hazardous waste. Do not allow reagent to reach surface water or sewage system.

Dispose of in accordance with official regulations.

Non-contaminated and completely empty packaging can be recycled.