Magnesium

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Magnesium

Reagent for quantitative In-vitro-determination of Magnesium in serum / plasma / urine of animals

MG 013

Order No. MG 013

Contents: 2 x 50 mL Colour reagent
Additionally needed:
Calcium-Magnesium-Standard Order No. CA ST
Diaglobal-Round Cuvettes, 40 pcs. Order No. LH 075

Method

Xylidyl blue^{1) 2)}

Sample material

Serum, heparinised plasma or urine
Do not use EDTA-plasma, haemolytical plasma or sera.
Add some drops of HCl to urine, dilute 1+4 with dest. Water.
Stability in serum: at +2°C to +25°C: 7 days

Reagent

Contents / concentrations of the colour reagent:
Xylidylblue 0.1 mmol/L; GEDTA 40 µmol/L; Potassium carbonate 70 mmol/L; Detergent 1.5%; Tris-buffer 200 mmol/L;
Sodium azide < 0.1%; Triton X-100 < 1%, pH 11,0

Safety information

The reagent contains sodium azide (< 0.1 %) and Triton X-100. Do not swallow and avoid contact with skin and mucous membranes.
A safety data sheet is available on our website.³⁾

Storage and shelf life

The colour reagent can be kept at a temperature between +2°C and +8°C until the expiry date indicated on the packaging.

Measurement conditions

Measurement device: VET Photometer Diaglobal

Meas. wavelength: 520 nm

Temperature: Room temperature

Measurement range

0.10 - 2.10 mmol/L (0.24 - 5.11 mg/dL)
Should values exceed this range, dilute sample 1+1 with physiological saline solution.
Multiply the result by 2.

Notes

1. Protect reagent from sunlight. Close the bottle after using.
2. The magnesium assay is very sensitive. We recommend the application of disposables.
3. For the "on-site" preparation of plasma we recommend the Diaglobal-Mini centrifuge (Order No. DZ 002)

Working instructions

Pipette into round cuvettes:			
	Blank	Standard	Analysis
Colour reagent	1500 µL	1500 µL	1500 µL
Standard	--	10 µL	--
Sample	--	--	10 µL
Mix well and measure after 5 minutes at room temperature standard and analysis against blank.			

- Select the <MG> test
- Insert cuvette with reagent into photometer (blank value)
- Remove blank cuvette after signal ton
- Insert cuvette with standard into photometer (standard-absorption is shown on the display)
- Remove cuvette
- Insert cuvette with sample (Analysis) into photometer
- Read the result

Tip

It is possible to prepare several vials for measurement at the same time by pre-pipetting 1500 µL of colour reagent into the appropriate number of Diaglobal round cuvettes.

Shelf life of the pipetted reagents in the vial: 5 hours

Calculation

Concentration of magnesium in serum/plasma/urine:
 $c \text{ [mmol/L]} = 0.8226 \times \text{Abs. analysis} / \text{Abs. standard}$
 $c \text{ [mg/dL]} = 2.0 \times \text{Abs. analysis} / \text{Abs. standard}$

Conversion in mg/dL: mmol/L x 2.4312 = mg/dL

Quality assurance

For quality assurance we recommend universal control sera from company Roche, www.roche.de:
PreciControl ClinChem Multi 1 / Multi 2 (4 x 5 mL)
Order-No.: 05 947 626 190 / 05 947 774 190

Reference values for cattle⁴⁾

	mmol/L	mg/dL
In serum / plasma	0.8 - 1.3	1.9 - 3.2

Summary⁴⁾

Magnesium is a component of many different enzymes which are involved in energy metabolism in the muscle cells. A reduced concentration of magnesium in the blood (hypomagnesaemia) leads to tetanic seizures, while an increase (hypermagnesaemia) leads to flaccid paralysis.

Evaluation⁴⁾

Hypermagnesaemia
- Parturient paresis
- Renal insufficiency
- Exsiccosis
- Latrogenic oversupply

Hypomagnesaemia

- Chronic bowel diseases (malabsorption)
- Low-magnesium diet
- Acute pancreatitis
- Cholestasis
- Renal insufficiency
- Hypercalcaemia
- Hypocalcaemic ketoacidosis
- Hyperthyroidism
- Tetany (cow, sheep)

Measurement principle^{1) 2)}

In alkaline solution with xylidyl blue, magnesium forms a colour complex which is measured at 520 nm. Calcium ions are masked with Glycol Ethylenediaminetetraacetic Acid (GEDTA).

The assessment is carried out using a standard.

$\text{Mg}^{2+} + \text{Xylidyl blue} \rightarrow \text{Magnesium colour complex}$

Performance parameters

Specificity / interferences

No interference caused by calcium ions.
Bilirubin up to 20 mg/dL and lipidaemia do not interfere.
Do not use haemolytic serums as erythrocytes have high magnesium concentrations.

Inaccuracy

The reproducibility was checked using human and control samples.

In series [n = 20]	Average [mmol/L]	Standard deviation [mmol/L]	VK [%]
Sample 1	0.88	0.02	2.3
Sample 2	1.15	0.02	2.1
Sample 3	1.45	0.02	1.6
From day to day [n = 20]	Average [mmol/L]	Standard-deviation [mmol/L]	VK [%]
Sample 1	0.67	0.04	5.8
Sample 2	0.93	0.05	5.4
Sample 3	1.38	0.05	3.9

Analytic sensitiveness

Lower detection limit: 0.08 mmol/L (0.2 mg/dL)

Comparison of methods

A comparison of the Diaglobal test MG 013 (y) and a commercially available test (x) resulted in the following correlation according to the linear regression:

$$y = 0.982x + 0.02$$
$$r = 0.993$$

n = 26, concentration range: 0.5 - 2.0 mmol/L

Information on disposal

Waste code number 180106:

Vials with reagent are considered hazardous waste. Do not allow reagent to reach surface water or sewage system.

Dispose of in accordance with official regulations.

Non-contaminated and completely empty packaging can be recycled.

Non-contaminated and completely empty packaging can be recycled.

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