**Best. Nr.** CK 121  
**Inhalt:** 20 Tests

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|                       |        |                      |
|-----------------------|--------|----------------------|
|                       | CK     |                      |
| Creatinphosphat + ADP | →      | Creatin + ATP        |
|                       | HK     |                      |
| Glucose + ATP         | →      | Glucose-6-P + ADP    |
|                       | G6P-DH |                      |
| Glucose-6-P + NADP    | →      | Gluconat-6-P + NADPH |

HK = Hexokinase, G6P-DH = Glucose-6-Phosphat-dehydrogenase, ADP = Adenosin-diphosphat, ATP = Adenosin-triphosphat

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**Safety information****Hazard statements**

H360D: May damage the unborn child

EUH208: Contains Imidazolidinyl urea. May produce an allergic reaction.

**Precautionary statements**

P201: Obtain special instructions before use.

P202: Do not handle until all safety precautions have been read and understood.

P263: Avoid contact during pregnancy and while nursing.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P308+P313: IF exposed or concerned: Get medical advice/attention.

P405: Store locked up.

A safety data sheet is available on request.<sup>3)</sup>

**Storage and shelf life**

The reagents have to be kept at a temperature between +2°C and +8°C until the expiry date indicated on the packaging.

Stability of ready-to-use reagent:

at +2°C to +8°C 48 h

at +20°C to +25°C 8 h

**Measurement conditions**

Measurement device: Vario Photometer II Diaglobal

Meas. wavelength: 365 nm

Temperature: 37°C

Additionally needed:

Mini centrifuge, Dry block heater (37°C)

**Measurement range**

Up to 1000 U/L

In case of exceeding this value, dilute the sample 1 + 4 with physiological saline solution, multiply the result by 5.

**Reference values**

Male < 170 U/L

Female < 145 U/L

**Preparation**

The dry block heater should be activated at least 20 min. before measurement.

**Working instructions**

|                                                           |          |
|-----------------------------------------------------------|----------|
| Pipette into round cuvettes:                              |          |
|                                                           | Analysis |
| Enzyme substrate solution                                 | 500 µL   |
| Pipette into ready-to-use reagent:                        |          |
| Serum / plasma                                            | 20 µL    |
| Mix thoroughly and put immediately into dry block heater. |          |

Up to six samples can be measured at the same time.

- Pipette 500µL enzyme substrate solution in cuvette
- with buffer solution (=ready-to-use reagent)
- Select <CK 121>
- Pipette 20µL serum / plasma in cuvette with ready-to-use reagent, mix well
- Insert cuvette into 37°C dry block heater and start measurement at the photometer
- Follow the operator guidance showing on display
- After finishing the measurements, read results in (U/L) with arrow keys
- Remove cuvette

**Calculation**

Activity (U/L) of creatine kinase in plasma will be calculated.

CK (U/L) = F x ΔA

The computation formula is stored in the Photometer.

**Quality assurance**

For quality assurance we recommend universal control

sera from company Roche, www.roche.de:

PreciControl ClinChem Multi 1 / Multi 2 (4 x 5 mL)

Order-No.: 05 947 626 190 / 05 947 774 190

Ref.: Roche / Hitachi analyzers, Method: IFCC liquid

**Advice**

For CK measurement in capillary blood, we recommend our test CK 321.

**Summary****Biochemical fundamentals**

Creatine Kinase has a key role in the catabolism of muscle cells and is responsible for reduction of ADP to ATP. Measurable total activity of CK in serum consists of activity of iso enzyme CK-MB, CK-MM and CK-BB. In blood plasma of a healthy person a very low concentration of CK is normally found. At muscle deterioration CK transgresses into blood and blood CK level increases. Since evacuation from cell occurs by lymph channel the concentration climax will be reached mostly after 6 - 8 hours.

**Training and competition control**

Increased CK level due to physical activity is an expression of overstraining or unaccustomed strain of individual muscle groups.

If the CK value increases to 900 U/L (15 µmol/Ls) during training, a decrease of training exercise is recommended or else overtraining is risked.<sup>4)</sup>

After extreme physical stress (marathon, triathlon) CK value can increase beyond 4000 U/L (70 µmol/Ls).<sup>4)</sup>

**Measurement principle**

With Diaglobal test CK 121 total CK is measured.

Determination is based on a 37°C standard method of IFCC (International Federation of Clinical Chemistry).<sup>5)</sup>

HK = Hexokinase, G6P-DH = Glucose-6-Phosphate-dehydrogenase, ADP = Adenosine-diphosphate, ATP = Adenosine-triphosphate

The increase of NADPH is directly proportional to the CK activity and is measured at 365 nm.

Note: Older technical literature often gives values that were determined with the help of the older 25°C method. Values, measured with the new IFCC method, average higher about a factor of 1.6.

**Performance parameters****Specificity / interferences<sup>6)</sup>**

Falsely increased values due to strong haemolysis (Hb >2.0 g/L), no significant interferences due to bilirubin (up to 15 mg/dL).

Strong lipemic samples can cause problems due to high absorptions. IN this case dilute the sample 1 + 4 with physiological saline solution.

**Inaccuracy**

The reproducibility was checked using human and control samples.

| In series<br>[n = 20]       | Average<br>[U/L] | Standard deviation<br>[U/L] | CV<br>[%]  |
|-----------------------------|------------------|-----------------------------|------------|
| Sample 1<br>Sample 2        | 156<br>502       | 4.84<br>9.54                | 3.1<br>1.9 |
| From day to day<br>[n = 20] | Average<br>[U/L] | Standard deviation<br>[U/L] | CV<br>[%]  |
| Sample 1<br>Sample 2        | 152<br>499       | 5.32<br>12.0                | 3.5<br>2.4 |

**Analytic sensitiveness**

Lower detection limit: 20 U/L (0.33 µkat/L)

**Comparison of methods**

A comparison of the Diaglobal test CK 121 (y) and another commercial available test (x) resulted in the following correlation according to the linear regression:

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