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
HK = Hexokinase, G6P-DH = Glucose-6-Phosphat-dehydrogenase,
 ADP = Adenosin-diphosphat, ATP = Adenosin-triphosphat

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Reagent for quantitative In-vitro-determination of CK in blood

CK 321

Order no. CK 321
Contents: 20 tests
20 pipette tips, 500µL

Method
IFCC 37°C^{1,2)} / DGKC modified

Sample material
Capillary or EDTA blood.
Set capillary blood immediately in reaction tube "R".
Stability of the sample in reaction tube "R":
6 hours (+15°C to +25°C)

Reagents
Contents / concentrations:
1. Starter reagent (caps in PE-bottle)
ADP 2 mmol/L, G6-PDH >2 kU/L
2. Buffer solution (pre-portioned in round cuvettes)
Imidazole buffer 100 mmol/L, pH = 6.00 (37°C), Glucose 20 mmol/L,
Magnesium acetate 10 mmol/L, EDTA 2 mmol/L, AMP 5 mmol/L,
NADP 2 mmol/L, N-Acetyl-cysteine 20 mmol/L, Diadenosin-penta-
phosphate 10 µmol/L, HK > 4.0 kU/L; Sodium azide < 0.1 %
3. Solution for dilution of blood (reaction tubes "R")
Creatine phosphate 30 mmol/L, Sodium chloride 3 g/L,
Sodium azide < 0.1 %

Safety information
Hazard statements
H360D: May damage the unborn child
EUH208: Contains Imidazolidinyl urea. May produce an allergic reaction.
Precautionary statements
P201: Obtain special instructions before use.
P202: Do not handle until all safety precautions have been read and understood.
P263: Avoid contact during pregnancy and while nursing.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P308+P313: IF exposed or concerned: Get medical advice/attention.
P405: Store locked up.

A safety data sheet is available on request.³⁾

Storage and shelf life
The reagent have to be kept in a dark place at a temperature between + 2°C to + 8°C until the expiry date indicated on the packaging.

Measurement conditions
Measurement device: Vario Photometer II Diaglobal
Meas. wavelength: 365 nm
Temperature: 37°C
Additionally required:
Mini centrifuge, Dry block heater (37°C)

Measurement range
Up to 1000 U/L
In case of exceeding these value, repeat determination with sample volume of 20 µL (instead of 60 µL). Multiply the result by 3

Reference values
Male < 170 U/L
Female < 145 U/L

Preparation
Reagents are ready to use.
The dry block heater should be activated at least 20 min. before measurement.

Working instructions
Up to six samples can be measured at the same time.

- Withdraw 60µL capillary blood from finger pad or earlobe with end-to-end capillary and insert into reaction tube "R"
- Transfer blood into solution by mixing strongly
- Insert reaction tube "R" into Minicentrifuge and centrifuge for 1 minute
- Pipette 500µL supernatant into round cuvettes (blanks)
- Select <CK 321>
- Screw caps from PE-bottle onto prepared round cuvettes and mix very strongly to solve the starter reagent
- Insert cuvettes into 37°C dry block heater and start measurement at the same time
- Follow the operator guidance showing on display
- Upon completion of measurement: Adjust HCT value with arrow keys and confirm with <ENTER> *)
- After finishing the measurements, read results in U/L with arrow keys
- Remove cuvette

**) HCT value of 40 % is set as standard value.
Please note that in case of measuring CK 321 with control serum the HCT-value has to be keyed in as 0%.*

Calculation
Activity (U/L) of Creatine kinase in plasma will be calculated. For exact measurement results the haematocrit values has to be considered.

$$CK (U/L) = F \times \Delta E / (1 - 0.01 \times HCT)$$

F = calibration constant, HCT = Haematocrit in %
The algorithm for calculating the results is stored in the Vario Photometer II.


Conversion: µkat/L x 60 = U/L

Quality assurance
For quality assurance we recommend universal control sera from company Roche, www.roche.de:
PreciControl ClinChem Multi 1 / Multi 2 (4 x 5 mL)
Order-No.: 05 947 626 190 / 05 947 774 190
Ref.: Roche / Hitachi analyzers, Method: IFCC liquid

Summary
Biochemical fundamentals
Creatine Kinase has a key role in the catabolism of muscle cells and is responsible for reduction of ADP to ATP. Measurable total activity of CK in serum consists of activity of iso enzyme CK-MB, CK-MM and CK-BB. In blood plasma of a healthy person a very low concentration of CK is normally found. At muscle deterioration CK transgresses into blood and blood CK level increases. Since evacuation from cell occur by lymph channel the concentration climax will be reached mostly after 6 - 8 hours.

Training and competition control
Increased CK level due to physical activity is an expression of overstraining or unaccustomed strain of individual muscle groups.
If the CK value increases to 900 U/L (15 µmol/Ls) during training, a decrease of training exercise is recommended or else overtraining is risked.⁴⁾
After extreme physical stress (marathon, triathlon) CK value can increase beyond 4000 U/L (70 µmol/Ls).⁴⁾

Measurement principle
With Diaglobal test CK 321 total CK is measured.
Determination is based on a 37°C standard method of IFCC⁵⁾ (International Federation of Clinical Chemistry).

HK = Hexokinase, G6P-DH = Glucose-6-Phosphate-dehydrogenase,
ADP = Adenosine-diphosphate, ATP = Adenosine-triphosphate

The increase of NADPH is directly proportional to the CK activity and is measured at 365 nm.

Note: Earlier technical literature often gives values that were determined with the help of the 25°C method. Values, measured with the new IFCC method, average higher about a factor of 1.6.

Performance parameters
Specificity / interferences⁶⁾
Falsely increased values due to strong haemolysate (Hb>2.0 g/L), no significant interferences due to bilirubin (up to 15 mg/dL).
Strong lipemic samples can cause problems due to high absorptions.

Inaccuracy
The reproducibility was checked using human and control samples.

In series [n = 20]	Average [U/L]	Standard deviation [U/L]	CV [%]
Sample 1	150	4.96	3.3
Sample 2	503	12.3	2.4
From day to day [n = 20]	Average [U/L]	Standard deviation [U/L]	CV [%]
Sample 1	148	5.48	3.7
Sample 2	508	14.2	2.8

Analytic sensitiveness
Lower detection limit: 20 U/L (0.33 µkat/L)

Comparison of methods
A comparison of the Diaglobal test CK 321 (y) and another commercial available test (x) resulted in the following correlation according to the linear regression:
 $y = 0.984x - 3.65$
 $r = 0.996$
 $n = 45$
Concentration range: 40 - 2600 U/L

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